



## Transformation and Production of Chemically Competent Cells

www.davids-bio.com

### Buffers needed for the Procedure

WASH BUFFER:	STORAGE BUFFER:	10 x SOC	LB MEDIUM
100 mM RbCl	10 mM RbCl	0.2 % (w/v) KCL	0.5 % (m/v) Yeast Extract
50 mM MnCl <sub>2</sub>	75 mM CaCl <sub>2</sub>	2.0 % (w/v) MgCl <sub>2</sub>	1 % (m/v) Tryptone
30 mM KOAc	10 mM MOPS	2.0 % (w/v) MgSO <sub>4</sub>	1 % (m/v) NaCl
10 mM CaCl <sub>2</sub>	15 % Glycerol	4.0 % (w/v) Glucose	(1.5 % (m/v) Agar for plates only)
15 % Glycerol		In <i>A. dest</i> (sterile filtered)	In <i>A. dest</i> (autoclave)

### Production of Chemically Competent Cells (D. Hanahan, J. Mol. Biol., 1983, 166, 557-580.)

- Create an Overnight Culture (LB medium, 37°C, 200 rpm) of your desired *E. coli* cell (i.e. DH5α)
- Inoculate 100 ml main culture 1 : 100 with the overnight culture (1 ml overnight culture in 100 ml LB)
- Cultivate the cells at 37°C and 200 rpm until the OD<sub>600</sub> = 0.3 – 0.4 is reached
- Chill the culture on ice for 10 min
- Centrifuge the cold cells at 4000 x g at 4°C and discard the supernatant
- Resuspend the pellet in 20 ml Wash Buffer and incubate the cells on ice for 15 min
- Centrifuge the cold cells at 4000 x g at 4°C and discard the supernatant
- Resuspend the cells in 4 ml Storage Buffer and incubate the cells on ice for 15 min
- Aliquot the cells in your desired volume (i. e. 50 µl each aliquot)
- Shock freeze the aliquots immediately with liquid nitrogen or a “Dry Ice – Water – NaCl” solution
- Store the cells at -80°C until use

### Chemical Transformation

- Thaw an aliquot of the chemical competent cells on ice for 5 min
- Add 1 µl plasmid (10 – 100 ng) to the cells and stir the tube with the pipet
- Incubate the cells on ice for 30 min
- Heat shock the cells exactly 30 seconds at 42°C (use a water bath)
- Chill the cells on ice for 2 min
- Add 450 µl 1 x SOC (dilute the 10x SOC 1:10 in LB)
- Incubate the cells at 37°C and 200 rpm for 1 h
- Plate the cells on an agar plate with the desired antibiotic (you may use different dilutions of your transformation approach to gain single colonies)